

Manipulating Small Volumes of Glycerol

~ WEAR GLOVES ~

I. Practicing with a P-20

1. Label two empty reaction tubes A and B with a permanent ink marker, if not already labeled.
2. Add the amounts of solutions I, II, III and V to tubes A and B as shown in the table below.

Solution Volume in Each Tube						
Reaction Tube	10% glycerol	25% glycerol	40% glycerol	60% glycerol	80% glycerol	Total in Tube
A	4 μ L	5 μ L	2 μ L			
B	6.5 μ L	2.5 μ L			2 μ L	

3. Spin tubes A and B in the microcentrifuge for 1-2 sec to pool the solutions.

> REFER TO PREVIOUS CENTRIFUGE INSTRUCTIONS

4. Add up the total volume of liquid in Tube A. AS A CHECK OF YOUR TECHNIQUE, set the micropipet to that volume and withdraw all of the liquid in tube A. The contents should *just* fill the tip -- no air space at the bottom of the tip; no leftover liquid in the tube. Discard liquid and tip into waste beaker
5. How much should be in tube B? μ L Check your technique by setting the pipettor to the correct volume and withdraw all (hopefully) the solution in tube B.

II. Practicing with a P-200

Reaction Tube	10% glycerol	25% glycerol	40% glycerol	60% glycerol	80% glycerol	Total in Tube
C	20 μ L		57 μ L	110 μ L		

1. Label an empty reaction tube C, if not already labeled.
2. Add the amounts of solutions I, II and IV to tube C as shown in the table below.
3. Spin tube C for 1-2 sec. (Be sure to balance your tube.)
4. Check the accuracy of your micropipeting technique with the P-200: Set the pipet to μ L (the predicted "supposed to be" volume) & withdraw the contents of tube C.

III. Practicing with a P-1000.

1. Label an empty reaction tube D if not already labeled
2. Add the amounts of solutions III and V to tube D as shown in the table below.

Reaction Tube	10% glycerol	25% glycerol	40% glycerol	60% glycerol	80% glycerol	Total in Tube
D	100 μ L		320 μ L		580 μ L	

- Spin tube D for 1-2 sec. (Be sure to balance your tube.)
- Check the accuracy of your micropipeting technique with the P-1000. Set the pipet to _____ μL and withdraw the contents of tube D.

IV. Practicing with all micropipets

- Label an empty reaction tube E, if not already labeled.
- Fill in the blanks in the following chart. Indicate on the line "P-_____" the appropriate micropipet to be used for each sample. Then write in the three boxes below each line, the numbers that should be dialed to give you the indicated volume.
- Add the volumes of solutions I-V indicated above into tube E, using the *appropriate micropipet*.

Reaction Tube E	10% glycerol	25% glycerol	40% glycerol	60% glycerol	80% glycerol
Volume:	15 μl	105 μl	2 μl	38 μl	350 μl
Micropipet: (P-20,P-200,P-1000)	P-_____	P-_____	P-_____	P-_____	P-_____
Micropipet setting					

Check the accuracy of your micropipeting technique with the P-_____. Set the pipet to _____ μL and withdraw the contents of tube E.